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Item	Item	Hours
1	質子治療系統開機前檢查 (Morning check before start up PTS)	5
2	迴旋加速器開機程序 (Start-up cyclotron)	3
3	迴旋加速器與射束傳輸系統各項元件之操作 (Operation of each equipment of cyclotron and beam transport system)	20
4	迴旋加速器射束輸出品質調控 (Beam tuning of cyclotron)	20
5	射束傳輸至治療室之操作 (Beam transportation to treatment room)	10
6	迴旋加速器之關機程序 (Shutdown cyclotron)	2
7	設備作業之每日備份 (Recordation of daily operation)	5
8	故障排除 (Trouble shooting)	30
9	離子源之例行維護 (Maintenance of ion source)	20
10	真空系統之例行維護 (Maintenance of vacuum pumps)	10
11	其他設備之例行維護 (Maintenance of other equipment)	10
12	消耗性零件與備用零件之管理 (Control of consumable and spare parts)	5
13	電力系統中斷後之關機與重開機程序 (Shutdown and start-up procedure for planned electric power interruption)	10
14	意外事故 (電力中斷、地震、火災等) 中斷治療後的復歸與檢查作業 [Check and recovery procedure after accident (power failure, earthquake, fire etc.)]	10
	總時數	160

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	校 驗 項 目	
校驗頻率	ICRU 59 (1998)	ICRU 78 (2007)
Daily	Lasers and radiographic systems	Aperture alignment; room lasers, interlocks; communication; patient-positioning system
	Beam range for each acceleration energy	Depth-dose and lateral profiles (range, entrance dose, uniformity of range modulation and Bragg-peak width, flatness, symmetry)
	Dose monitor calibration for specific sets of standard conditions	Dose monitor calibration, check of MU value under standard condition
	Individual patient treatment calibration and range check	Individual patient treatment calibration and range checks
	Uniformity of range modulation	
	Backup dose monitor function	

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	校 驗 項 目	
校驗頻率	ICRU 59 (1998)	ICRU 78 (2007)
Weekly	Dose delivered for a random selection of sample patient treatments	Patient-positioning and imaging systems
	Treatment field lateral uniformity	Beam-line apparatus
	Room interlock checks	Respiratory-gating equipment
		Dose delivered to randomly selected patients (comparison of planned dose distributions to those measured in a water phantom)

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	校 驗 項 目	
校驗頻率	ICRU 59 (1998)	ICRU 78 (2007)
Annual	Patient alignment systems	x-ray patient positioning and alignment systems
	Gantry isocenter location and uncertainty	CT Hounsfield number calibration
	Patient support system motions and readouts	Comprehensive tests of therapy equipment:
	Depth dose distribution for monoenergetic and modulated beams	1. monitor chambers, timers, beam-delivery termination and control interlocks, stray radiation exposure to patients, gantry.
	Depth dose distribution for monoenergetic and modulated beams	2. isocenter, depth-dose and lateral profiles, baseline data for daily QA checks
	Lateral uniformity and penumbra measurements for a range of conditions	
	Dependence of dose monitor calibration on treatment parameters	
	Base line data for daily quality assurance checks	
	Calibration intercomparison with other proton facilities	

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ICRU 59 (1998)	
校驗頻率	校 驗 項 目
New treatment portals	Dose monitor calibration
	Consistency of planned treatment range with beam energy and absorber thicknesses
	Coincidence of the physical aperture shape with the planned shape
	Qualitative consistency of the range compensator shape with treatment plan isothickness lines
	Comparison of range compensator physical thicknesses measured at selected points with planned thickness
	Agreement of dose monitor calibration with predictions from a computational model

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ICRU 78, Additional quality-assurance procedures for scanning beam-delivery systems	
校驗頻率	校 驗 項 目
Daily	Dose rate and monitor ratios for the pencil beam
	Performance of the beam-position monitors
	Depth-dose curve of a pencil beam in a water phantom
	Calibration of the primary dose monitor
Weekly	Qualitative three-dimensional check of the outline and range of the dose distribution for one patient' s irradiation field in a water phantom
Half-yearly	Calibration of the primary dose monitor and the phase space of the beam tunes
Annual	Check of the beam characteristics:
	calibration of the whole dosimetry system, performance of the scanning system in terms of dose linearity and dose-rate dependence

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Daily QA for Scattering and Uniform Scanning Beams

Test	Tolerance/Comments
Review operator's cyclotron and gantry startup checklists	Operator checks machine operating parameters daily.
Safety interlocks, indicator lights, neutron detector, A/V systems	
Output constancy check for scattering and uniform scanning reference fields	2%. Output measurements in plastic phantom.
Dose rate measurement	30% tolerance.
kV imaging and laser accuracy	1 mm. Orthogonal x-ray crosshair and laser agreement.
Range verifier measurements of scattering and uniform scanning fields	1 mm.
Range modulation wheel control signal timing	0.2 ms. Review of system reported timing values compared to baseline data.
Scan field size length and width	3 mm. Compare system reported values to expected values for uniform scanning reference field.

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Weekly QA for Scattering and Uniform Scanning Beams

Test	Tolerance/Comments
Review daily QA results	Review for system performance trends.
MLIC calibration	MLIC calibration performed prior to use for PDD measurements.
Output constancy check for two scattering and two uniform scanning reference fields	2%. Output measurements using MLIC. Fields measured include intermediate and large range fields.
Range measurements of above fields	1.5 mm. Measurements performed using MLIC.
SOBP width measurements of above fields	2 mm. Measurements performed using MLIC.
SOBP uniformity measurements of above fields	2%. Measurements performed using MLIC.

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Monthly QA for Scattering and Uniform Scanning Beams

Test	Tolerance/Comments
Review daily/weekly QA results	Review for system performance trends.
Output constancy check for two scattering and two uniform scanning reference fields plus random patient field	2%. Output measurements performed in water. Fields measured include intermediate and large range fields. Patient field measured to support MU calculation accuracy.
Range measurements of above fields	1.5 mm. Measurements performed using in-water scans.
SOBP width measurements of above fields	2 mm. Measurements performed using in-water scans.
SOBP uniformity measurements of above fields	2%. Measurements performed using in-water scans.
Lateral dose profile measurements of two scattering and two uniform scanning reference fields	3% flatness; 1.5% symmetry; 2-mm field size. Measurements performed using MatriXX ion chamber array.
Range verifier reading of pristine peak ranges with and without first scatterers in beam	1 mm. First scatterer integrity tested.
X-ray/light field/proton beam field agreement	1 mm. Double-exposure film measurements. Rotating through snouts.
Patient-positioning system/gantry accuracy	1 mm/1°. Tests include isocentric table rotation.
Digital imaging positioning system	1 mm in calculated vs. set shifts. Tests performed for two gantry angles.

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Annual QA for Scattering and Uniform Scanning Beams

Test	Tolerance/Comments
Review daily/weekly/monthly QA results	Review for system performance trends.
Dosimetry system intercomparisons	All dosimetry systems used for daily, weekly, and monthly QA tests are compared to ADCL-calibrated dosimetry system and calibration factors determined.
Output calibration	1.5%. Scattering and uniform scanning reference field outputs calibrated using ADCL-calibrated dosimetry system.
MU chamber linearity	0.5%
MU chamber dependence on dose rate	1% from 1 to 3 Gy/min
Dose rate accuracy	25% for scattering and uniform scanning reference fields.
Output factors	2%. Performed for each option.
Output dependence on gantry angles	0.5%. Output of reference fields measured at prime gantry angles.
Depth-dose scans	1 mm for range, 3 mm for SOBP width, 2% for SOBP uniformity. Measurements performed for each option. Off axis scans performed for 1 field.
Lateral dose profile scans	2% for flatness and symmetry. Measurements performed for each option.
Lateral dose profile dependence on gantry angle	0.5% variations.
Gantry accuracy and isocentricity	0.5° and -mm diameter.
Patient-positioning system translation, rotation accuracy, and isocentricity	0.5 mm, 0.2°, and 1-mm diameter.
Snout position accuracy and alignment	5 mm in snout position accuracy; 0.5 mm in alignment in both in-plane and cross-plane directions.
Imaging system accuracy	Within 1 mm for all snouts in both in-plane and cross-plane directions.
Safety interlocks and radiation monitors	A number of safety interlocks and in-room neutron monitors are tested.

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1. Particle Therapy Statistics in 2013. International Journal of Particle Therapy, Martin Jermann, MSc, 2014.
2. ICRU REPORT 78, Prescribing, Recording, and Reporting Proton-Beam Therapy. 2007.
3. ICRU (1998). Clinical Proton Dosimetry - Part I Beam Production Beam Delivery and Measurement of Absorbed Dose. ICRU REPORT 59. BETHESDA, MARYLAND, U.S.A, International Commission on Radiation Units and Measurements.